

Clino Guard – easy installation, reliable protection

On a daily basis, nursing staff and policy makers in nursing homes and hospitals are confronted with the two terms cost pressure and staff shortage. This is contrasted with directly or indirectly involved people who expect the best care, a lot of comfort and most importantly security.

This is not an easy job - but with intelligent, combined systems those procedural challenges can be solved efficiently.



Honeywell Life Safety wants to meet these market requirements with Clino Guard from Ackermann by Honeywell and offers a RFID based system to warn the nursing staff if a patient is leaving the defined security area.

Complete security for the patient

Clino Guard unifies security with free moving space, especially for patients that suffer from dementia. If a patient who has a transponder is leaving the defined area, the system starts special activities: for example it sends an alarm to the responsible nursing staff, actuates the observation camera activates the functions that keep the doors closed. Furthermore the system knows automatically if the patient is going out with a staff member and in this case it doesn't activate the alarm.

The sending/receiving units can get installed everywhere to allow unremarkable protection: in the door area, in indoor rooms or outside the building. The range of the system is between 3,5 and 6 meters, you can increase it if you want with more antenna and contact plates. We provide reliable recognition and data transfer by using different rest and working frequencies and three-dimensional acting antenna.



Besides the three-dimensional antenna the patient transponder are also equipped with observation of the battery and are only able to open with a special key.

Your advantages with Clino Guard

- Reliable security for patient and staff
- Unremarkable because of the invisible security areas
- you do not need much effort to service the system
- more life quality and free moving space
- robust, discreet transponder for patient
- transponder for the staff that allows to stop the alarm
- easy installation of the system

Complete security for the staff

Besides the security actions for the patients, we also provide security to the staff, which is possible because of the secure logging of all calling information (place, time, person and confirmation), the connection to the nurse call system and the access monitoring of entrances, elevators, etc.

The fitting system for everybody

Clino Guard is available in two variants. The simple Clino Guard system provides the observation of smaller buildings and area with small effort. For bigger projects and connections with other systems we suggest Clino Guard plus

Inform yourself on our website www.hls-austria.com or call us directly: 01/6006030-0